TeSeE™ WESTERN BLOT

32 TESTS

355 1169

**REAGENTS FOR *in vitro* CONFIRMATION OF SUSPECTED
TSE POSITIVE SAMPLES**

Fitness for purpose validated and certified by OIE
Registration number: 20090105

BIO-RAD

CONTENTS

- 1 - GENERAL INFORMATION
- 2 - ASSAY PRINCIPLE
- 3 - COMPOSITION OF THE KIT
- 4 - SAMPLES
- 5 - ASSAY PROCEDURE WITH MINI BLOT™ GEL
 - 5.1 Additional reagents and material required
 - 5.2 Preparation of reagents
 - 5.3 Sample purification
 - 5.4 Electrophoresis
 - 5.5 Protein transfer
 - 5.6 Immunoblotting
- 6 - ASSAY PROCEDURE WITH CRITERION™ XT GEL
 - 6.1 Additional reagents and material required
 - 6.2 Preparation of reagents
 - 6.3 Sample purification
 - 6.4 Electrophoresis
 - 6.5 Protein transfer
 - 6.6 Immunoblotting
- 7 - INTERPRETATION OF RESULTS
- 8 - PRECAUTIONS
- 9 - HYGIENE AND SAFETY INSTRUCTIONS
- 10 - REFERENCES

1 - GENERAL INFORMATION

Transmissible Spongiform Encephalopathies (TSE's) were first reported in the eighteenth century in sheep (Scrapie) and more recently in cervids such as deer and elk (Chronic Wasting disease, CWD) and cattle (Bovine Spongiform Encephalopathy, BSE). Humans are also susceptible to certain forms of TSE such as Kuru, Creutzfeldt-Jakob Disease (CJD) or Gerstmann-Sträussler-Scheinker Syndrome (GSS). The emergence of new variant Creutzfeldt-Jakob Disease (vCJD) in the human population has been strongly linked to the dietary intake of BSE-infected meat or meat products. One of the main characteristics of TSEs is a progressive accumulation in the central nervous system of an abnormal isoform of natural or cellular prion protein (PrP^C), termed PrP^{res}. This disease specific PrP^{res} is characterised by an increased resistance to proteases. The TeSeE™ WESTERN BLOT permits qualitative identification of PrP^{res} after proteolytic treatment which results in a reduced molecular weight fragment due to 'N' terminus truncation.

Active/passive surveillance programs have been conducted worldwide to detect BSE, scrapie or CWD in infected animals. Those programs have resulted in the identification of increased numbers of positive cases at the screening laboratories. Those positive samples (suspected animals) are then systematically confirmed as "TSE-infected" by the demonstration of typical spongiform changes with histopathology, or with the detection of abnormal PrP by Immunohistochemistry (IHC), or of Scrapie Associated Fibrils (SAFs) by electron microscopy. These above confirmation techniques require technical expertise for the interpretation of the results and are time consuming and expensive. Western blot technique can also be considered as an alternative method for confirmation of the TSE suspected samples.

The TeSeE™ WESTERN BLOT is using the same assay principle than the Bio-Rad rapid assays (TeSeE™ SAP, TeSeE™ sheep/goat) that includes the preliminary purification and concentration of the PrP^{res}, associated to a highly sensitive immunoblotting. Then, it can be used efficiently for confirmation of any TSE suspected samples and for typing of TSE strains in sheep.

2 - ASSAY PRINCIPLE

The TeSeE™ WESTERN BLOT assay allows the detection of PrP^{res} in nervous tissues (bovine, ovine, caprine, cervids, ...) or peripheral tissues (cervids) collected from infected animals.

The assay procedure begins with the digestion of cellular PrP protein (PrP^c), followed by purification and concentration of disease specific PrP^{res}. Detection of PrP^{res} is carried out by electrophoretic migration then immunoblotting using a monoclonal antibody highly specific for PrP^{res}.

The assay procedure includes the following steps:

- Sample homogenization,
- Digestion of PrP^c with proteinase K,
- Purification and concentration of PrP^{res},
- Electrophoresis and transfer onto a membrane,
- Immunoblotting.

3 - COMPOSITION OF THE KIT

| Labelling | Type of reagents | Presentation | Storage |
|----------------|---|---------------------|---------------|
| Grinding Tubes | Grinding tubes containing ceramic beads in a buffer solution ⁽¹⁾ | 1 bag (35 tubes) | +2°C to +25°C |
| A | Denaturing solution Ready to use | 1 vial (20 ml) | +2°C to +25°C |
| B | Clarifying solution Colouring: bromophenol blue Ready to use | 1 vial (20 ml) | +2°C to +25°C |
| PK | Proteinase K Colouring: phenol red | 1 vial (0.5 ml) | +2°C to +8°C |
| Ab I | Primary antibody ⁽¹⁾ : anti-PrP monoclonal antibody (10x) | 1 vial (8 ml) | +2°C to +8°C |
| Ab II | Secondary antibody ⁽¹⁾ : Sheep Anti-Mouse IgG-(H+L)-HRP (10x) | 1 vial (10 ml) | +2°C to +8°C |
| Bl | Blocking solution ⁽¹⁾ (10x) | 1 vial (10 ml) | +2°C to +8°C |

⁽¹⁾ These reagents contain 0.1 % of ProClin® 300 (preservative).

4 - SAMPLES

TeSeE™ WESTERN BLOT assay is suitable for the detection of TSEs in cattle (Bovine Spongiform Encephalopathy, BSE), in ovine and caprine (BSE and scrapie), and in cervids (Chronic Wasting Disease, CWD).

This test can be processed directly from the same sample homogenate (grinding tube) prepared for Bio-Rad rapid testing (TeSeE™ SAP, TeSeE™ sheep/goat).

Bovine: purification of PrP^{res} is performed on samples from Central Nervous System (CNS). Since distribution of PrP^{res} is heterogeneous in central nervous system, obex area from brainstem must be preferably sampled for optimal detection.

Sampling syringe (Ref.: 355 1175) allows easy and rapid sampling of obex area in a secure way. Please refer to sampling protocol provided with sample syringes for detailed instructions on good sampling procedure.

Small ruminants: purification of PrP^{res} is performed on samples from Central Nervous System (CNS).

Samples are cut and weighed individually.

Cervids: purification of PrP^{res} is performed on samples from Central Nervous System (CNS) or peripheral tissues (lymphoid nodes).

Samples are cut and weighed individually.

5 - ASSAY PROCEDURE WITH MINI BLOT™ GEL

5.1 - ADDITIONAL REAGENTS AND MATERIAL REQUIRED

5.1.1 - REAGENTS AND DISPOSABLES

Graduated pipettes (5, 10, 25 ml), conical tubes (50 ml), 2 ml polypropylene micro-test tubes with caps.

PARAFIL™M Sealing films.

Sample purification

| | | |
|-----------------------|-------|----------------------------|
| Laemmli sample buffer | 30 ml | Bio-Rad, cat. Nr. 161-0737 |
| 2-Mercaptoethanol | 25 ml | Bio-Rad, cat. Nr. 161-0710 |
| SDS | 100 g | Bio-Rad, cat. Nr. 161-0301 |
| Calibration syringes | 200 | Bio-Rad, cat. Nr. 355-1174 |

Gel electrophoresis

| | | |
|---|--------|-----------------------------|
| Acrylamide 40% 29:1 | 500 ml | Bio-Rad, cat. Nr. 161-0146 |
| 0.5 M Tris-HCl pH 6.8 | 1 L | Bio-Rad, cat. Nr. 161-0799 |
| 1.5 M Tris-HCl pH 8.8 | 1 L | Bio-Rad, cat. Nr. 161-0798 |
| Bromophenol blue | 10 g | Bio-Rad, cat. Nr. 161-0404 |
| Sucrose | 1 kg | Bio-Rad, cat. Nr. 161-0720 |
| Ammonium persulfate | 10 g | Bio-Rad, cat. Nr. 161-0700 |
| TEMED | 5 ml | Bio-Rad, cat. Nr. 161-0800 |
| Tris/Glycine/SDS (running buffer) (10x) | 1 L | Bio-Rad, cat. Nr. 161-0732 |
| Kaleidoscope prestained standard | 500 µl | Bio-Rad, cat. Nr. 161-0324 |
| MagicMark™ XP Western Standard (Molecular weight standard) | 250 µl | Invitrogen, cat. Nr. LC5602 |

Immunoblotting

| | | |
|---|-----------|----------------------------|
| Ethanol (Normapur) | 1L | VWR, cat. Nr. 20821-296 |
| Tris/CAPS (transfer buffer) (10x) | 1L | Bio-Rad, cat. Nr. 161-0778 |
| Filter paper (transfer paper for Mini Blot™ handcast gels) | 50 sheets | Bio-Rad, cat. Nr. 170-3932 |
| PVDF membrane (0.2 µm) | 10 sheets | Bio-Rad, cat. Nr. 162-0175 |
| Tween® 20 | 100 ml | Bio-Rad, cat. Nr. 170-6531 |
| PBS (washing buffer) (10x) | 1 L | Bio-Rad, cat. Nr. 161-0780 |

| | | |
|--------------------------------|------------|--------------------------------------|
| ECL (substrate for conjugate) | 125 ml | Amersham, cat. Nr. RPN2109 |
| ECL Hyperfilms (18 x 24 cm) | 25 films | Amersham, cat. Nr. RPN2103K |
| Development folders | 30 folders | Applied Biosystems, cat Nr. T2258 |
| Kodak developing solution LX24 | to 20 L | VWR or Kodak |
| Kodak fixative solution AL4 | to 20 L | VWR or Kodak |

5.1.2 - EQUIPMENT

Adjustable pipettes (10, 40, 200, 1000 µl),
 Graduated cylinder (1 L and 2 L), plastic forceps, trays, Vortex®.
 Exposure cassette and red light for film development.

Sample purification

| | |
|----------------------------------|----------------------------|
| TeSeE™ PRECESS 48™ | Bio-Rad, cat. Nr. 359-0200 |
| TeSeE™ PRECESS 24™ | Bio-Rad, cat. Nr. 359-1070 |
| Or Ribolyser® | Bio-Rad, cat. Nr. 358-9158 |
| Heating block | Bio-Rad, cat. Nr. 358-9046 |
| Heating block adaptor - 20 tubes | Bio-Rad, cat. Nr. 358-9199 |
| Centrifuge - 220/240 V | Bio-Rad, cat. Nr. 358-9190 |
| Drum rotor | Bio-Rad, cat. Nr. 358-9189 |
| Rotor adaptors - (x6) | Bio-Rad, cat. Nr. 358-9191 |

Gel electrophoresis

| | |
|--|----------------------------|
| Mini-PROTEAN® Tetra Cell, electrophoresis module | Bio-Rad, cat. Nr. 165-8002 |
| 5 spacers plates | Bio-Rad, cat. Nr. 165-3312 |
| PowerPac HC power supply: 100/120 V - 220/240 V | Bio-Rad, cat. Nr. 164-5052 |

Transfer

| | |
|------------------|----------------------------|
| Trans-Blot® Cell | Bio-Rad, cat. Nr. 170-3946 |
|------------------|----------------------------|

Immunoblotting

| | |
|----------------------------------|----------------------------|
| Western Processor base | Bio-Rad, cat. Nr. 359-0098 |
| Western Processor Mini Blot™ kit | Bio-Rad, cat. Nr. 170-3988 |

5.2 - PREPARATION OF REAGENTS

5.2.1 - SAMPLE PURIFICATION

• Proteinase K

Solution of proteinase K diluted in reagent A:

- ▶ 1 ml Reagent A
- ▶ 20 µl Proteinase K

Mix well by inverting until you obtain a homogeneous solution. After reconstitution, diluted proteinase K is stable 10 hours at room temperature (+18°C to +30°C).

• Laemmli solution

Solution of SDS + 2-Mercaptoethanol + Laemmli sample buffer:

- ▶ 0.6 g SDS
- ▶ 1.5 ml 2-Mercaptoethanol

Mix by inverting.

- ▶ 28.5 ml Laemmli sample buffer

Solution is aliquoted into 4 ml aliquots and stored at -20°C. Thawed aliquots can be re-frozen.

Note: It is recommended to prepare Laemmli solution one hour before use allowing SDS to be completely dissolved.

5.2.2 - ELECTROPHORESIS

• Hand cast discontinuous acrylamide gel

The gel must be 1.5 mm thickness.

Using the Mini Blot™ casting module, the resolving gel (13.5% acrylamide, pH 8.8) is cast first, once the resolving gel is polymerized the stacking gel is added (3% acrylamide, pH 6.8).

Resolving gel (1 gel)

- ▶ 2.8 ml Acrylamide 40%, 29:1
- ▶ 1.7 ml 1.5 M Tris-HCl buffer, pH 8.8 / SDS (1)
- ▶ 1.3 ml 50% sucrose solution (2)
- ▶ 2.5 ml distilled water

Mix by inverting.

- ▶ 43 µl 10% Ammonium persulfate (3)
- ▶ 9 µl TEMED

Pour 7 ml of the gel solution into the plates and retain the residual solution as a control of polymerization. Gently overlay to the top with 1 ml of 0.3 M Tris-HCl pH 8.8 / SDS buffer (4) so that the gel surface doesn't dry out. Let the gel polymerize for 15-20 minutes at room temperature (+18°C to +30°C). Check the residual solution is polymerized. Invert the plate assembly to eliminate excess of buffer.

Stacking gel (1 gel)

- ▶ 4 ml 3% Acrylamide solution (7)
- ▶ 28 µl 10% Ammonium persulfate (3)
- ▶ 6 µl TEMED

Mix by inverting.

Gently pour the stacking gel onto the resolving gel and retain the residual solution as a control of polymerization. Position the comb, taking care not to trap any bubble in the well positions.

Let the gel polymerize for 5-10 minutes at room temperature (+18°C to +30°C). Check the residual solution is polymerized.

(1) Solution of 1.5 M Tris-HCl buffer, pH 8.8 / SDS

- ▶ 0.2 g SDS
- ▶ 50 ml 1.5 M Tris-HCl buffer pH 8.8

Solution can be stored at +2°C to +8°C for 2 weeks.

(2) Solution of 50% Sucrose

- ▶ 25 g Sucrose
- ▶ to 50 ml Distilled water

Sucrose solution can be stored at +2°C to +8°C for 2 weeks.

(3) Solution of 10% Ammonium persulfate

- ▶ 5 g Ammonium persulfate
- ▶ to 50 ml Distilled water

Ammonium persulfate solution is aliquoted and stored at -20°C. Thawed solution can be stored at +2°C to +8°C for 2 weeks.

(4) Solution of 0.3 M Tris-HCl buffer, pH 8.8 / SDS

- ▶ 40 ml Distilled water
- ▶ 10 ml 1.5 M Tris-HCl buffer pH 8.8 / SDS

Solution can be stored at +2°C to +8°C for 2 weeks.

(5) Solution of 0.5 M Tris-HCl buffer, pH 6.8 / SDS

- ▶ 0.2 g SDS
- ▶ 50 ml 0.5 M Tris-HCl buffer pH 6.8

Solution can be stored at +2°C to +8°C for 2 weeks.

(6) Solution of 1% Bromophenol Blue

- ▶ 0.5 g Bromophenol Blue
- ▶ 50 ml Distilled water

Bromophenol Blue solution can be stored at room temperature (+18°C to +30°C) for 6 months.

(7) Solution of 3% acrylamide

- ▶ 3.8 ml Acrylamide 40%, 29:1
- ▶ 10 ml 0.5 M Tris-HCl buffer pH 6.8 / SDS (5)
- ▶ 6 ml Sucrose 50% (2)
- ▶ 500 µl Bromophenol Blue 1% (6)
- ▶ to 50 ml Distilled water

Solution can be stored at +2°C to +8°C for 2 weeks.

• Kaleidoscope prestained standard

The kaleidoscope prestained standard is prepared during the sample denaturation before loading on the acrylamide gel.

Prepare a 1/12 dilution in Laemmli solution (for example 10 µl of the kaleidoscope prestained standard + 110 µl of Laemmli solution).

Please refer to the kaleidoscope prestained standard insert for storage conditions.

• MagicMark™ XP Western Standard

The MagicMark™ XP molecular weight is prepared during the sample denaturation before loading on the acrylamide gel.

Prepare a 1/12 dilution in Laemmli solution (for example 10 µl of MagicMark™ XP + 110 µl of Laemmli solution).

Please refer to MagicMark™ XP insert for storage conditions.

• Mini Blot™ migration buffer

Solution of Tris-Glycine-SDS (1x).

Prepare a 1/10 dilution. **1 L of diluted buffer is required for 1 tank:**

- ▶ 900 ml Distilled water
- ▶ 100 ml Tris-Glycine-SDS buffer (10x)

Homogenize. Solution can not be stored.

5.2.3 - PROTEIN TRANSFER

• Transfer buffer

Solution of Tris/CAPS-Ethanol 15%. **2.5 L is required for 1 transfer tank.**

- ▶ 750 ml Distilled water
- ▶ 150 ml Pure ethanol
- ▶ 100 ml Tris/CAPS (10x)

Homogenize. Solution can not be stored.

5.2.4 - IMMUNOBLOTTING

• Wash solution 1

Solution of PBS (1x) + 0.1% Tween® 20. **Approximately 500 ml is required for the complete process of 1 membrane.**

- ▶ 900 ml Distilled water
- ▶ 100 ml PBS (10x)
- ▶ 1 ml Tween® 20

Thoroughly homogenize. Solution can be stored at +2°C to +8°C, overnight.

• Wash solution 2

Solution of PBS (1x). **Approximately 100 ml is required for the complete process of 1 membrane.**

- ▶ 900 ml Distilled water
- ▶ 100 ml PBS (10x)

Solution can be stored at room temperature (+18°C to +30°C) overnight.

• Blocking solution

During the transfer step, dilute the blocking solution (Bl) 1/10 in Wash solution 1. **20 ml of diluted blocking solution (1x) is required for 1 membrane.**

- ▶ 18 ml Wash solution 1
- ▶ 2 ml Blocking solution (10x)

Homogenize by tube inverting.

• Diluted primary antibody

Just prior to use, dilute the primary antibody 1/10 in Wash solution 1. **15 ml of diluted primary antibody is required for 1 membrane.**

- ▶ 13.5 ml Wash solution 1
- ▶ 1.5 ml Primary antibody (10x)

Homogenize by tube inverting.

- **Diluted secondary antibody (conjugate)**

Just before use, dilute the secondary antibody 1/10 in Wash solution 1.

- **20 ml of diluted conjugate is required for 1 membrane.**

- ▶ 18 ml Wash solution 1
- ▶ 2 ml Secondary antibody (10x)

Homogenize by tube inverting.

- **ECL**

Substrate (ECL) must be prepared just before use. **1 ml of substrate is required for 1 membrane.**

- ▶ 0.5 ml Reagent 1
- ▶ 0.5 ml Reagent 2

Homogenize the solution.

- **Development solution**

- ▶ 800 ml Distilled water
- ▶ 200 ml Development product

Solution can be stored at room temperature (+18°C to +30°C), in a darkroom for 15 days maximum.

- **Fixative solution**

- ▶ 800 ml Distilled water
- ▶ 200 ml Fixative product

Solution can be stored at room temperature (+18°C to +30°C), in a darkroom for 15 days maximum.

5.3 - SAMPLE PURIFICATION

TeSeE™ WESTERN BLOT can be processed directly from the same sample homogenate (grinding tube) prepared for Bio-Rad rapid tests (TeSeE™ SAP, TeSeE™ sheep/goat).

Sampling

For peripheral tissues (lymph nodes) insert one medium bead (Ref.: 355 1171) in the grinding tube, before to add the sample.

Take a mass of 350 mg ± 40 mg of nervous tissue (preferably in the obex area) or 200 mg ± 20 mg of peripheral tissue.

Deposit the sample in grinding tube, close firmly and proceed to the grinding step in the homogenizer (Ribolyser®, TeSeE™ PRECESS 48™ or TeSeE™ PRECESS 24™-system).

Sample grinding

Place the tubes in the crown of the homogenizer.

Perform one agitation cycle with the following instrument parameters.

| | Ribolyser® | | TeSeE™ PRECESS 48™ or TeSeE™ PRECESS 24™ | |
|-------------|-----------------|-----------------------|---|--------------------|
| | Nervous tissues | Peripheral tissues | Nervous tissues | Peripheral tissues |
| Time (sec.) | 45 | 2 x 45 ⁽¹⁾ | – | – |
| Speed | 6.5 | 6.5 | – | – |
| Program | – | – | Program 1 | Program 2 |

When grinding is insufficient, another 1 or 2 agitation cycles can be performed⁽²⁾.

⁽¹⁾ ⁽²⁾ Wait a 5 minutes pause between the 2 agitation cycles.

Sample calibration

Remove the grinding tubes from the homogenizer, resuspend the homogenate by inverting before opening the tubes and aspirate 500 µl with the calibration syringe taking care to immerse the needle below the level of ceramic beads to avoid sampling tissue fragments.

Transfer each 500 µl sample into a 2 ml Eppendorf micro test-tube.

Note: At this stage, both grinding tubes after homogenisation and micro test-tubes after sample calibration can be stored, closed:

- At room temperature (+18°C to +30°C) for 15 hours.
- At +2°C to +8°C for 72 hours.
- At -20°C for 1 year. Frozen samples must be thawed at room temperature (+18°C to +30°C).

Samples can be submitted to a maximum of 3 freezing/thawing cycles. Samples must always be homogenized by inverting before use.

Proteinase K treatment

Distribute 500 µl of reconstituted proteinase K solution (see paragraph 5.2.1) into each micro test-tube.

Homogenize the sealed tubes by inverting (10 times) and incubate at $37^{\circ}\text{C} \pm 2^{\circ}\text{C}$ in a heating block incubator for 10 minutes.

Precipitation of PrP^{res} with reagent B

Remove the tubes from the incubator. Open them and distribute 500 µl of reagent B into each tube. Homogenize by inverting until a homogeneous colour is obtained.

Concentration of the PrP^{res} by centrifugation

Centrifuge the tubes for 7 minutes at 15000 g at 20°C .

Sample clarifying

Discard the supernatant by inverting over a waste container. Then dry the tubes by inverting onto absorbent paper for 5 minutes.

Distribute 100 µl of the Laemmli solution (see paragraph 5.2.1) into each micro test-tube.

Incubate for 5 minutes at room temperature ($+18^{\circ}\text{C}$ to $+30^{\circ}\text{C}$).

Completely resolubilize the pellet by aspiration/dispensing with a pipette.

Incubate for 5 minutes at $100^{\circ}\text{C} \pm 5^{\circ}\text{C}$ in a heating block incubator.

Remove the tubes from the incubator, homogenize by vortexing.

Centrifuge the tubes for 15 minutes at 15000 g at 20°C .

Transfer the supernatant to a new micro test-tube. Discard the tube containing the pellet.

At this stage, the supernatant can be stored frozen at -20°C for 24 hours; the samples must be thawed at room temperature ($+18^{\circ}\text{C}$ to $+30^{\circ}\text{C}$) prior to use.

5.4 - ELECTROPHORESIS

The TeSeE™ WESTERN BLOT assay can be used for both confirmation of TSE suspected samples or for strain typing in sheep.

The following procedure is applicable for confirmation of TSE suspected samples.

Please contact your Bio-rad representative for instruction protocol in case of strain typing application.

Gel preparation

Place the acrylamide gels (see paragraph 5.2.2) in the migration tank. Pour the migration buffer (see paragraph 5.2.2) into the electrophoresis tank on each side of the gels, up to the top of the wells. Carefully remove the combs and rinse each well with migration buffer, using a pipette.

Sample loading

Heat the samples for 4 minutes at $100\text{ }^{\circ}\text{C} \pm 5\text{ }^{\circ}\text{C}$ just prior to loading 15 μl /well. Load 15 μl of the diluted kaleidoscope prestained standard and 15 μl of the diluted MagicMark™ XP (see paragraph 5.2.2).

Note: In case several gels are run at the same time, make sure that you stagger the loading of controls into different lanes for easy identification.

Differential migration of the samples

Run the gel at room temperature ($+18\text{ }^{\circ}\text{C}$ to $+30\text{ }^{\circ}\text{C}$) for 90 minutes at 150 V. The blue line must be out of the gel.

5.5 - PROTEIN TRANSFER

The transfer buffer must be prepared before the end of the sample migration (see paragraph 5.2.3).

Protein transfer preparation

Cut the membrane to the gel dimensions. Always use forceps when handling the membrane.

Immerse the membrane in pure ethanol for 15 seconds, rinse in distilled water for 5 minutes, then for 10 minutes in the transfer buffer.

Carefully remove the gel from the glass plates and let it equilibrate for 10 minutes in the transfer buffer.

Gel sandwich preparation

Soak filter paper and fibre pads in the transfer buffer.

Open the transfer cassette, with transparent side on the left. Respectively place on the transparent side a fibre pad, a filter paper, the membrane* and the gel*. Complete with a filter paper then a fibre pad and close the cassette.

Immerse it in the transfer tank previously filled to the indicated limit with transfer buffer.

*Remove any air bubbles which may have formed.

Note: In case several membranes are processed at the same time, label each membrane in the corner.

Transfer onto the PVDF membrane

Agitate during the transfer by using a magnetic stirring bar and run for 60 minutes at 115 V.

5.6 - IMMUNOBLOTTING

- a) Upon completion of the protein transfer, open the blotting assembly and remove the membrane for development. Quickly immerse the membrane in Wash solution 2 (see paragraph 5.2.4), then place it in ethanol for 10 seconds before rinsing for 5 minutes in distilled water.

Note: At this step, the membrane can be stored overnight in distilled water at +2°C to +8°C.

Let the membrane adjust to room temperature (+18°C to +30°C) before to start the immunoblotting.

- b) Eliminate distilled water and incubate the membrane for 30 minutes in blocking solution (see paragraph 5.2.4). Incubate under medium agitation.

20 ml is sufficient for 1 membrane.

Note: from this step until the step g), the Bio-Rad Western Processor can be used for agitation and washing steps (refer to instruction manual for settings).

- c) Eliminate the blocking solution and incubate the membrane in diluted **primary antibody** (see paragraph 5.2.4) for 30 minutes at room temperature (+18°C to +30°C) under medium agitation.

15 ml of diluted primary antibody is required for 1 membrane.

- d) Eliminate the primary antibody solution and using Wash solution 1, briefly rinse the membrane, then wash twice for respectively 5 and 10 minutes, under fast agitation.

50 ml of Wash solution 1 is required for each cycle and for 1 membrane.

- e) Eliminate Wash solution 1 and incubate the membrane for 20 minutes in diluted **secondary antibody** (see paragraph 5.2.4) at room temperature (+18°C to +30°C) under medium agitation.
20 ml of diluted secondary antibody is required for 1 membrane.
- f) Eliminate the secondary antibody solution and using Wash solution 1, briefly rinse, then wash for respectively 5, 10 and 10 minutes under fast agitation.
50 ml of Wash solution 1 is required for each cycle and for 1 membrane.
- g) Place the membrane in 50 ml of Wash solution 2 under slow agitation.
- h) Drain the membrane on absorbent paper without blotting and place it in the plastic folder.
- i) Add the ECL reagent (see paragraph 5.2.4). Eliminate the excess of reagent and air bubbles with absorbent paper. Place into the exposure cassette.
- j) In a darkroom, cover the folder with a film and expose for 15 minutes. Film can be exposed longer or shorter time for optimal signal.
- k) Immerse the film in developing solution for 45 seconds (see paragraph 5.2.4). Rinse in distilled water. Immerse the film in fixative solution until the film becomes completely transparent.
- l) Wash with distilled water and let the film dry.

6 - ASSAY PROCEDURE WITH CRITERION™ XT GEL

6.1 - ADDITIONAL REAGENTS AND MATERIAL REQUIRED

6.1.1 - REAGENTS AND DISPOSABLES

Graduated pipettes (5, 10, 25 ml), conical tubes (50 ml), 2 ml polypropylene micro-test tubes with caps.

PARAFILM®M Sealing films.

Sample purification

| | | |
|-----------------------|-------|----------------------------|
| Laemmli sample buffer | 30 ml | Bio-Rad, cat. Nr. 161-0737 |
| 2-Mercaptoethanol | 25 ml | Bio-Rad, cat. Nr. 161-0710 |
| SDS | 100 g | Bio-Rad, cat. Nr. 161-0301 |
| Calibration syringes | 200 | Bio-Rad, cat. Nr. 355-1174 |

Gel electrophoresis

| | | |
|--|------------------|-----------------------------|
| Criterion™ XT 12 % Bis-Tris | 1 gel - 18 wells | Bio-Rad, cat. Nr. 345-0118 |
| XT-MOPS (running buffer) (20 x) | 500 ml | Bio-Rad, cat. Nr. 161-0788 |
| Kaleidoscope prestained standard | 500 µl | Bio-Rad, cat. Nr. 161-0324 |
| MagicMark™ XP Western Standard (Molecular weight standard) | 250 µl | Invitrogen, cat. Nr. LC5602 |

Immunoblotting

| | | |
|--|------------|------------------------------------|
| Ethanol (Normapur) | 1L | VWR, cat. Nr. 20821-296 |
| Tris/CAPS (transfer buffer) (10x) | 1L | Bio-Rad, cat. Nr. 161-0778 |
| Filter paper (transfer paper for Criterion™ XT precast gels) | 50 sheets | Bio-Rad, cat. Nr. 170-4085 |
| PVDF membrane (0.2 µm) | 10 sheets | Bio-Rad, cat. Nr. 162-0175 |
| Tween® 20 | 100 ml | Bio-Rad, cat. Nr. 170-6531 |
| PBS (washing buffer) (10x) | 1 L | Bio-Rad, cat. Nr. 161-0780 |
| ECL (substrate for conjugate) | 125 ml | Amersham, cat. Nr. RPN2109 |
| ECL Hyperfilms (18 x 24 cm) | 25 films | Amersham, cat. Nr. RPN2103K |
| Development folders | 30 folders | Applied Biosystems, cat. Nr. T2258 |
| Kodak developing solution LX24 | to 20 L | VWR or Kodak |
| Kodak fixative solution AL4 | to 20 L | VWR or Kodak |

6.1.2 - EQUIPMENT

Adjustable pipettes (10, 40, 200, 1000 µl).

Graduated cylinder (1L and 2L).

Plastic forceps, Trays, Vortex®.

Exposure cassette and red light for film development.

Sample purification

| | |
|----------------------------------|----------------------------|
| TeSeE™ PRECESS 48™, | Bio-Rad, cat. Nr. 359-0200 |
| TeSeE™ PRECESS 24™, | Bio-Rad, cat. Nr. 359-1070 |
| Or Ribolyser® | Bio-Rad, cat. Nr. 358-9158 |
| Heating block | Bio-Rad, cat. Nr. 358-9046 |
| Heating block adaptor - 20 tubes | Bio-Rad, cat. Nr. 358-9199 |
| Centrifuge - 220/240 V | Bio-Rad, cat. Nr. 358-9190 |
| Drum rotor | Bio-Rad, cat. Nr. 358-9189 |
| Rotor adaptors - (x6) | Bio-Rad, cat. Nr. 358-9191 |

Gel electrophoresis

| | |
|--|----------------------------|
| Criterion™ XT Cell | Bio-Rad, cat. Nr. 165-6001 |
| PowerPac HC power supply: 100/120 V - 220/240 V | Bio-Rad, cat. Nr. 164-5052 |

Transfer

| | |
|-----------------------|----------------------------|
| Criterion™ XT blotter | Bio-Rad, cat. Nr. 170-4070 |
|-----------------------|----------------------------|

Immunoblotting

| | |
|-------------------------------------|----------------------------|
| Western Processor base | Bio-Rad, cat. Nr. 359-0098 |
| Western Processor Criterion™ XT kit | Bio-Rad, cat. Nr. 170-3985 |

6.2 - PREPARATION OF REAGENTS

6.2.1 - SAMPLE PURIFICATION

- **Proteinase K**

Solution of proteinase K diluted in reagent A:

- ▶ 1 ml Reagent A
- ▶ 20 μ l Proteinase K

Mix well by inverting until you obtain a homogeneous solution. After reconstitution, diluted proteinase K is stable 10 hours at room temperature (+18°C to +30°C).

- **Laemmli solution**

Solution of SDS + 2-Mercaptoethanol + Laemmli sample buffer:

- ▶ 0.6 g SDS
- ▶ 1.5 ml 2-Mercaptoethanol

Mix by inverting.

- ▶ 28.5 ml Laemmli sample buffer

Solution is aliquoted into 4 ml aliquots and stored at -20°C. Thawed aliquots can be re-frozen.

Note: It is recommended to prepare Laemmli solution one hour before use allowing SDS to be correctly dissolved.

6.2.2 - ELECTROPHORESIS

- **Kaleidoscope prestained standard**

The kaleidoscope prestained standard is prepared during the sample denaturation before loading on the acrylamide gel.

Prepare a 1/12 dilution in Laemmli solution, for example 10 μ l of the kaleidoscope prestained standard + 110 μ l of Laemmli solution.

Please refer to the kaleidoscope prestained standard insert for storage conditions.

- **MagicMark™ XP Western Standard**

The MagicMark™ XP molecular weight is prepared during the sample denaturation before loading on the acrylamide gel.

Prepare a 1/12 dilution in Laemmli solution, for example 10 μ l of MagicMark™ XP + 110 μ l of Laemmli solution.

Please refer to MagicMark™ XP insert for storage conditions.

• **Criterion™ XT migration buffer**

Solution of MOPS (1x).

Prepare a 1/20 dilution. **1 L of diluted buffer is required for 1 tank:**

- ▶ 950 ml Distilled water
- ▶ 50 ml MOPS buffer (20x)

Homogenize. Solution can not be stored.

6.2.3 - PROTEIN TRANSFER

• **Transfer buffer**

Solution of Tris/CAPS-Ethanol 15%. **Approximately 2 L is required for 1 migration tank.**

- ▶ 750 ml Distilled water
- ▶ 150 ml Pure ethanol
- ▶ 100 ml Tris/CAPS (10x)

Homogenize. Solution can not be stored.

6.2.4 - IMMUNOBLOTTING

• **Wash solution 1**

Solution of PBS (1x) + 0.1% Tween® 20. **Approximately 1 L is required for the complete process of 1 membrane.**

- ▶ 900 ml Distilled water
- ▶ 100 ml PBS (10x)
- ▶ 1 ml Tween® 20

Thoroughly homogenize. Solution is stored at +2°C to +8°C, overnight.

• **Wash solution 2**

Solution of PBS (1x). **Approximately 200 ml is required for the complete process of 1 membrane.**

- ▶ 900 ml Distilled water
- ▶ 100 ml PBS (10x)

Solution is stored at room temperature (+18°C to +30°C) overnight.

• **Blocking solution**

During the transfer step, dilute the blocking solution (Bl) 1/10 in Wash solution 1. **40 ml of diluted blocking solution is required for 1 membrane.**

- ▶ 36 ml Wash solution 1
- ▶ 4 ml Blocking solution (10x)

Homogenize by tube inverting.

- **Diluted primary antibody**

Just prior to use, dilute the primary antibody 1/10 in Wash solution 1.

30 ml of diluted antibody is required for 1 membrane.

- ▶ 27 ml Wash solution 1
- ▶ 3 ml Primary antibody (10x)

Homogenize by inverting.

- **Diluted secondary antibody (conjugate)**

Just prior to use, dilute the secondary antibody 1/10 in Wash solution 1.

40 ml of diluted conjugate is required for 1 membrane.

- ▶ 36 ml Wash solution 1
- ▶ 4 ml Secondary antibody (10x)

Homogenize by tube inverting.

- **ECL**

Substrate (ECL) must be prepared just prior to use. **2 ml of substrate is required for 1 membrane.**

- ▶ 1 ml Reagent 1
- ▶ 1 ml Reagent 2

Homogenize.

- **Development solution**

- ▶ 800 ml Distilled water
- ▶ 200 ml Development product

Solution is stored at room temperature (+18°C to +30°C), in a darkroom for 15 days maximum.

- **Fixative solution**

- ▶ 800 ml Distilled water
- ▶ 200 ml Fixative product

Solution is stored at room temperature (+18°C to +30°C), in a darkroom for 15 days maximum.

6.3 - SAMPLE PURIFICATION

TeSeE™ WESTERN BLOT can be processed directly from the same sample homogenate (grinding tube) prepared for Bio-Rad rapid tests (TeSeE™SAP, TeSeE™ sheep/goat).

Sampling

For peripheral tissues (lymph nodes) insert one medium bead (Ref.: 355 1171) in the grinding tube, before to add the sample.

Take a mass of 350 mg ± 40 mg of nervous tissue (preferably in the obex area) or 200 mg ± 20 mg of peripheral tissue.

Deposit the sample in grinding tube, close firmly and proceed to the grinding step in the homogenizer (Ribolyser®, TeSeE™ PRECESS 48™ or TeSeE™ PRECESS 24™ - system).

Sample grinding

Place the tubes in the crown of the homogenizer.

Perform one agitation cycle with the following instrument parameters.

| | Ribolyser® | | TeSeE™ PRECESS 48™ or TeSeE™ PRECESS 24™ | |
|-------------|-----------------|-----------------------|---|--------------------|
| | Nervous tissues | Peripheral tissues | Nervous tissues | Peripheral tissues |
| Time (sec.) | 45 | 2 x 45 ⁽¹⁾ | – | – |
| Speed | 6.5 | 6.5 | – | – |
| Program | – | – | Program 1 | Program 2 |

When grinding is insufficient, another 1 or 2 agitation cycles can be performed⁽²⁾.

^{(1) (2)} Wait a 5 minutes pause between the 2 agitation cycles.

Sample calibration

Remove the grinding tubes from the homogenizer, resuspend the homogenate by inverting before opening the tubes and aspirate 500 µl with the calibration syringe taking care to immerse the needle below the level of ceramic beads to avoid sampling tissue fragments.

Transfert each 500 µl sample into 2 ml Eppendorf micro test-tube.

Note: at this stage, both grinding tubes after homogenisation and micro test-tubes after sample calibration can be stored, closed:

- At room temperature (+18°C to +30°C) for 15 hours.
- At +2°C to +8°C for 72 hours.
- At -20°C for 1 year. Frozen samples must be thawed at room temperature (+18°C to +30°C).

Samples can be submitted to a maximum of 3 freezing/thawing cycles. Samples must always be homogenized by inverting before use.

Proteinase K Treatment

Distribute 500 µl of reconstituted proteinase K solution (see paragraph 6.2.1) into each micro test-tube.

Homogenize the sealed tubes by inverting (10 times) and incubate at $37^{\circ}\text{C} \pm 2^{\circ}\text{C}$ in a heating block incubator for 10 minutes.

Precipitation of PrP^{res} with reagent B

Remove the tubes from the incubator. Open and distribute 500 µl of reagent B into each tube. Homogenize by inverting until a homogeneous colour is obtained.

Concentration of the PrP^{res} by centrifugation

Centrifuge the tubes for 7 minutes at 15000 g at 20°C .

Sample clarifying

Discard the supernatant by inverting over a waste container. Then dry the tubes by inverting onto absorbent paper for 5 minutes.

Distribute 100 µl of the Laemmli solution (see paragraph 6.2.1) into each micro test-tube.

Incubate for 5 minutes at room temperature ($+18^{\circ}\text{C}$ to $+30^{\circ}\text{C}$).

Completely resolubilise the pellet by aspiration/dispensing with a pipette.

Incubate for 5 minutes at $100^{\circ}\text{C} \pm 5^{\circ}\text{C}$ in a heating block incubator.

Remove the tubes from the incubator, homogenize by vortexing.

Centrifuge the tubes for 15 minutes at 15000 g at 20°C .

Transfer the supernatant to a new micro test-tube. Discard the tube containing the pellet.

At this stage, the supernatant can be stored frozen at -20°C for 24 hours; the samples must be thawed at room temperature ($+18^{\circ}\text{C}$ to $+30^{\circ}\text{C}$) prior to use.

6.4 - ELECTROPHORESIS

The TeSeE™ WESTERN BLOT assay can be used for both confirmation of TSE suspected samples or for strain typing in sheep.

The following procedure is applicable for confirmation of TSE suspected samples.

Please contact your Bio-rad representative for instruction protocol in case of strain typing application.

Gel preparation

Remove the plastic band on the bottom of the plastic plate and place the acrylamide gels (see paragraph 6.2.2) in the migration tank. Pour the migration buffer (see paragraph 6.2.2) on each side of the gel up to the top of the wells and into the electrophoresis tank. Carefully remove the combs and rinse each well with migration buffer, using a pipette.

Sample loading

Heat the samples for 4 minutes at $100\text{ }^{\circ}\text{C} \pm 5^{\circ}\text{C}$ just prior to loading 15 μl /well. Load 15 μl of the diluted kaleidoscope prestained standard and 15 μl of the diluted MagicMark™ XP (see paragraph 6.2.2).

Note: In case several gels are run at the same time, make sure that you stagger the loading of controls into different lanes for easy identification.

Differential migration of the samples

Run the gel at room temperature ($+18^{\circ}\text{C}$ to $+30^{\circ}\text{C}$) for 50 minutes at 200 V.

6.5 - PROTEIN TRANSFER

The transfer buffer must be prepared before the end of the sample migration (see paragraph 6.2.3).

Protein transfer preparation

Cut the membrane to the gel dimensions. Always use forceps when handling the membrane.

Immerse the membrane in pure ethanol for 15 seconds, rinse in distilled water for 5 minutes, then for 10 minutes in the transfer buffer.

Carefully remove the gel from the plastic plates and let it equilibrate for 10 minutes in the transfer buffer.

Gel sandwich preparation

Soak filter paper and fibre pads in the transfer buffer.

Open the transfer cassette, with red side on the left. Respectively place on the red side a fibre pad, a filter paper, the membrane* and the gel*.

Complete with a filter paper then a fibre pad and close the cassette.

Immerse it in the transfer tank, previously filled to the indicated limit with transfer buffer. A frozen ice pack is added prior to fill the tank.

*Remove any air bubbles which may have formed.

Note: In case several membranes are processed at the same time, label each membrane in the corner.

Transfer onto the PVDF membrane

Agitate during the transfer by using a magnetic stirring bar and run for 60 minutes at 115 V.

6.6 - IMMUNOBLOTTING

- a) Upon completion of the protein transfer, open the blotting assembly and remove the membrane for development. Quickly immerse the membrane in Wash solution 2 (see paragraph 6.2.4), then place it in ethanol for 10 seconds before rinsing for 5 minutes in distilled water.

Note: At this step, the membrane can be stored overnight in distilled water at +2°C to +8°C.

Let the membrane adjust to room temperature (+18°C to +30°C) before to start the immunoblotting.

- b) Eliminate distilled water and incubate the membrane for 30 minutes in blocking solution (see paragraph 6.2.4). Incubate under medium agitation.

40 ml is required for 1 membrane.

Note: from this step until step g), the Bio-Rad Western Processor can be used for agitation and washing steps (refer to instruction manual for settings).

- c) Eliminate the blocking solution and incubate the membrane in diluted **primary antibody** (see paragraph 6.2.4) for 30 minutes at room temperature (+18°C to +30°C) under medium agitation.

30 ml of diluted primary antibody is required for 1 membrane.

- d) Eliminate the primary antibody solution and using Wash solution 1, briefly rinse the membrane, then wash twice for respectively 5 and 10 minutes, under fast agitation.


100 ml of Wash solution 1 is required for each cycle and for 1 membrane.

- e) Eliminate the Wash solution 1 and incubate the membrane for 20 minutes in diluted **secondary antibody** (see paragraph 6.2.4) at room temperature (+18°C to +30°C) under medium agitation.
40 ml of diluted secondary antibody is required for 1 membrane.
- f) Eliminate the secondary antibody solution and using Wash solution 1, briefly rinse for respectively 5, 10 and 10 minutes under fast agitation.
100 ml of Wash solution 1 is required for each cycle and for 1 membrane.
- g) Place the membrane in 100 ml of Wash solution 2 under slow agitation.
- h) Drain the membrane on absorbent paper without blotting and place it in the plastic folder.
- i) Add the ECL reagent (see paragraph 6.2.4). Eliminate the excess of reagent and air bubbles with absorbent paper. Place into the exposure cassette.
- j) In a darkroom, cover the folder with a film and expose for 15 minutes. Film can be exposed longer or shorter time for optimal signal.
- k) Immerse the film in developing solution for 45 seconds (see paragraph 6.2.4). Rinse in distilled water. Immerse the film in fixative solution until the film becomes completely transparent.
- l) Wash with distilled water and let the film dry.

7 - INTERPRETATION OF RESULTS

Figure 1 shows the expected band patterns for TSE negative samples, TSE positive samples in various animal species and molecular weight controls (positions 3 and 9).

Figure 1

Negative samples (positions 1 and 10) were treated with proteinase K. They do not show any signal, since PrP^c was fully digested.

Positive samples were also all treated with proteinase K.

BSE positive bovine sample (position 2), **classical scrapie positive sample** (position 6), and **CWD positive sample** (position 4) show a typical 3 band pattern, demonstrating digestion of PrP^c and transformation of the disease specific prion protein into a proteinase resistant core fragment with reduced molecular mass following removal of the N-terminus part of the protein. The two higher bands correspond to mono and di-glycosylated forms (27-30 kDa) while the lower band corresponds to the non glycosylated form.

Ovine sample experimentally infected with BSE (position 5) is presenting a higher signal on the di-glycosylated than on the mono-glycosylated band. Nevertheless, this typical glycoprofile can not be considered as a sufficient proof of the infection of the animal with BSE. According to the Community Reference Laboratory (CRL) recommendations, a differentiation assay must be performed on this type of sample to conclude between scrapie and BSE. Please contact Bio-Rad for more information on the Bio-Rad Discriminatory test.

Atypical scrapie (e.g. Nor98) affected ovine sample (position 8) is presenting an atypical glycoprofile. A lower band is visible at approximately 12 kDa, while other upper bands are not located at the same positions compare to "typical" scrapie cases. Signal is also stronger on the lowest band than on the upper band.

Gel reading must be considered cautiously since a strong positive sample detected with the TeSeE™ WESTERN BLOT may hide the nearest negative or low positive sample.

Limits of the test:

A negative result means that the test sample does not contain detectable PrP^{res} by TeSeE™ WESTERN BLOT assay. However, as very low levels of PrP^{res} may not be detected, such a result does not exclude the possibility of infection.

Any sample with a negative result according to the TeSeE™ WESTERN BLOT assay interpretation criteria, following a positive rapid test result, should be tested with one of the other OIE certified confirmatory methods, Immunohistochemistry (IHC) or SAF-Immunoblot.

Any sample with a reproducible positive result according to the test interpretation criteria must be verified in accordance with current legal regulation.

8 - PRECAUTIONS

The quality of the data obtained depends on compliance with the following good laboratory practices:

- Reagents must be stored at the appropriate temperature (refer to supplier's indications).
- Do not use reagents whose shelf-life has expired.
- Do not use reconstituted proteinase K after 10 hours storage at room temperature (+18°C to +30°C).
- Do not mix or combine reagents derived from different batches of the TeSe™ WESTERN BLOT during the same manipulation, with the exception of grinding tubes, reagent A, reagent B and proteinase K.
- Allow the reagents and buffers to adjust to room temperature (+18°C to +30°C) for 30 minutes before use.
- Thoroughly reconstitute reagents, avoiding any contamination.
- Do not perform the test in the presence of reactive vapors (acids, alkalines, aldehydes) or dust, which could alter the enzymatic activity of the conjugate.
- The enzymatic reaction is very sensitive to all metals or metallic ions. Consequently, no metallic element must be in contact with the conjugate.
- Only use polypropylene tubes.
- Use clean glassware, rinsed in distilled water, or preferably disposable material.
- Use a new pipette tip for each sample.
- When starting electrophoresis and transfer, check that the 2 electrodes are in contact with buffer.
- All the rinsing times must be respected to avoid any excess background noise during final staining with ECL reagent.

9 - HYGIENE AND SAFETY INSTRUCTIONS

Generally, hygiene conditions, biosafety measures and good laboratory practices must be in agreement with the recommendations of national regulatory authorities.

- All reagents of the kit are intended for use in “*in vitro*” diagnosis.
- Wear disposable gloves when handling reagents and samples and wash your hands thoroughly after handling them.
- Do not pipette by mouth.
- Use polypropylene containers to avoid broken glass.
- All the materials directly in contact with the samples and the wash solutions must be considered as contaminated.
- Avoid splashing samples or solutions containing samples.
- Contaminated surfaces must be cleaned with 20 000 ppm sodium hypochlorite solution (bleach). When the contaminating liquid is an acid, contaminated surfaces must be first neutralized with sodium hydroxide before using bleach. Surfaces must be rinsed with distilled water, dried with ethanol and wiped with absorbent paper. The material used for cleaning must be discarded in a specific container for contaminated waste.
- Samples, material and contaminated products must be eliminated after decontamination:
 - either by soaking in 1 M sodium hydroxide (final concentration) for at least 1 hour at room temperature (+18°C to +30°C),
 - or by soaking in 20 000 ppm sodium hypochlorite solution for at least 1 hour at room temperature (+18°C to +30°C),
 - or by autoclaving at 134°C minimum for at least 18 minutes, under 3 bars of pressure.

Note: never autoclave solutions containing bleach or reagent B.

- All operations involved in Transmissible Spongiform Encephalopathy (TSE) screening tests are subject to regulations and must be performed in an isolated, limited and controlled access laboratory devoted exclusively to this activity. A laboratory coat, overshoes, gloves, mask with visor or simple mask with safety glasses are required to ensure the operator's safety.
- Operators must receive specific training concerning the risks related to TSEs agents or prions and the validated modes of decontamination for

unconventional agents. Biosafety measures must be in agreement with recommendations of regular authorities of the country.

- Neutralize and/or autoclave all wash solutions or wash wastes or any liquid containing biological samples prior to their elimination.
- Reagent B is a dangerous substance classified as nocive (> 25% alcohol) according to European legislation.
- Reagents containing 0.1% ProClin® 300 are classified as irritating preparations according to European legislation.

Xn
(Alcohol > 25%)
(0.1% ProClin® 300)

R : 10-22-37/38-41-43-67 Flammable. Harmful if swallowed. Irritating to respiratory system and skin. Risk of serious damage to eyes. May cause sensitisation by skin contact. Inhalation of vapour may cause drowsiness and dizziness.

S : 7/9-13-26-28-37/39-46 Keep container tightly closed and in a well ventilated place. Keep away from food, drink and animal feed. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin wash immediately with plenty of water. Wear suitable protecting clothings, gloves and eye/face protection. If swallowed, seek medical advice immediately and show this container or label.

10 - REFERENCES

1. S.B. PRUSINER (1991)
Molecular biology of prion diseases - *Science* 252:1515-1522.
2. J.B. KATZ, J.G. PEDERSEN, A.L. JENNY and W.D. TAYLOR (1992)
Assessment of Western immunoblotting for the confirmatory diagnosis of ovine scrapie and bovine spongiform encephalopathy (BSE) - *Journal of Veterinary Diagnostic Investigations*, 4, 447-449.
3. G.A.H. WELLS, Y.I. SPENCER and M. HARITANI (1994)
Configuration and topographic distribution of PrP in the central nervous system in bovine spongiform encephalopathy: an immunohistochemical study. In: *Slow Infections of the Central Nervous System*, J. BJORNSSON, R.I. CARP, A. LOVE and M. WISNIEWSKI - Eds *The New York Academy of Sciences*, pp 350-352.
4. J.P. DESLYS, E. COMOY, S. HAWKINS, S. SIMON, H. SCHIMMEL, G. WELLS, J. GRASSI and J. MOYNAGH (2001)
Screening slaughtered cattle for BSE - *Nature*: 409; 476-477.
5. S. BENESTAD, P. SARRADIN, B. THU, J. SCHÖNHEIT, M. TRANULIS and B. BRATBERG (2003)
Cases of scrapie with unusual features in Norway and designation of new type, Nor98. *Veterinary Record* 153, 202-208.
6. A. BUSCHMANN, G. LÜHKEN, J. SCHULTZ, G. ERHARDT and M.-H. GROSCHUP (2004)
Neuronal accumulation of abnormal prion protein in sheep carrying a scrapie-resistant genotype (PrP^{ARR/ARR}). *Journal of General Virology* 85, 2727-2733.
7. L. ORGE, A.GALO, C.MACHADO, C. LIMA, C. OCHOA, J. SILVA, M. RAMOS and J.-P. SIMAS (2004)
Identification of putative atypical scrapie in sheep in Portugal. *Journal of General Virology* 85, 3487-3491.
8. A. BUSCHMANN, A.-G. BIACABE, U. ZIEGLER, A. BENCSIK, J.-Y. MADEC, G. ERHARDT, G. LÜHKEN, T. BARON and M.-H. GROSCHUP (2004)
Atypical scrapie cases in Germany and France are identified by discrepant reaction patterns in BSE rapid tests. *Journal of Virological Methods* 117, 27-36.

9. H. DE BOSSCHERE, S. ROELS, S. BENESTAD, E. VANOPDENBOSCH (2004). Scrapie case similar to Nor98 diagnosed in Belgium via active surveillance. *Veterinary Record* 155, 707-708.
10. H. ONNASCH, H. M. GUNN, B.J. BRADSHAW, S. BENESTAD, H.F. BASSETT (2004). Two Irish cases of scrapie resembling Nor98. *Veterinary Record* 155, 636-637.
11. S. BENESTAD, P. SARRADIN, J.-M. BILHEUDE, J. GRASSI, H. LAUDE, O. ANDREOLETTI, T. MOLDAL and B. BRATBERG. Are there gold standard methods for the diagnosis of TSE ? Nor98 scrapie cases: atypical cases and their challenging diagnosis. Second International Chronic Wasting Disease Symposium, Madison - Wisconsin, USA.
12. BIACABE, A-G., LAPLANCHE, J.L., RYDER, S. AND BARON, T. (2004). Distinct molecular phenotypes in bovine prion diseases, *EMBO Reports* 5, 110-114.
13. CASALONE, C., ZANUSSO, G., ACUTIS, P., FERRARI, S., CAPUCCI, I., TAGLIAVINI, F., MONACO, S. and CARAMELLI, M. (2004). Identification of a second bovine amyloidotic spongiform encephalopathy: Molecular similarities with sporadic Creutzfeldt-Jakob disease. *PNAS* 101, 3065-70.
14. BUSCHMANN A; GRETZSCHEL A; BIACABE A-G; SCHIEBEL K; CORONA C; HOFFMANN C; EIDEN M; BARON T; CASALONE C; GROSCHUP M H (2006). Atypical BSE in Germany-proof of transmissibility and biochemical characterization. *Veterinary microbiology* 2006;117(2-4):103-16.
15. ARSAC, J.-N., ANDREOLETTI, O., BILHEUDE, J.-M., LACROUX, C., BENESTAD, S. L., AND BAARON, T. (2007). Similar biochemical signatures and prion protein genotypes in atypical scrapie and Nor98 cases, France and Norway. *Emerging Infectious Diseases* 13(1), 58-65.
16. E. URO-COSTE, H CASSARD, S. SIMON, S. LUGAN, J.-M. BILHEUDE, A. PERRET-LIAUDET, J.-W. IRONSIDE, S. HAIK, C. BASSET-LEOBON, C. LACROUX, K. PEOCH, N. STREICHENBERGER, J. LANGEVELD, M.-W. HEAD, J. GRASSI, J.-J. HAUW, F. SCHELCHER, M.-B. DELISLE, O. ANDREOLETTI (2008) Beyond Pr^{Pres} Type 1/Type 2 Dichotomy in Creutzfeldt-Jakob disease. *Plos Pathogens* volume 4, issue 2, e1000029.

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